

Bock Manual lugano



Tubular Steel Products
Interior Furnishings
Therapy and Rehabilitation Beds

bock

Assembly and
Operating Instructions


Dear Customer,

In deciding to buy a therapy and rehabilitation bed from Bock you have opted for a care product that has a long service life and delivers first class functionality at the highest safety level. Our electrically adjustable care beds guarantee optimum comfort when lying, and support professional care activities. The focus is on people who need care, encouragement and protection. We have created the basic requirements for this with our care products. We urge you to prevent potential malfunction and risk of accidents by complying strictly with the safety and operating instructions and carrying out the necessary maintenance.

Yours sincerely



Klaus Bock



Content

4	1. General information
	1.1 <i>Practical – no packaging</i>
	1.2 <i>First impression – visual inspection</i>
6	2. Cleaning, care and disinfection
	2.1 <i>Cleaning and care</i>
	2.2 <i>Disinfection</i>
	2.3 <i>Avoiding danger</i>
9	3. General description of function
	3.1 <i>Structure and function</i>
	3.2 <i>How to assemble fast and easily</i>
	3.3 <i>Caution: Risk of injury</i>
17	4. Electric components
	4.1 <i>Drive unit</i>
	4.2 <i>The Control Box</i>
	4.3 <i>Height adjustment drive unit</i>
	4.4 <i>Lockable hand control, single fault safety</i>
	4.5 <i>Caution: Electrical drive unit</i>
21	5. Drive units
	5.1 <i>OKIMAT drive systems</i>
	5.2 <i>Mains isolation</i>
24	6. Accessories
	6.1 <i>Special sizes</i>
	6.2 <i>Assembling bed extensions</i>
	6.3 <i>Assembling accessories</i>
	6.4 <i>Mattresses</i>
	6.5 <i>Special mattress ripomed</i>
29	7. Setting up and operating – lugano
	7.1 <i>Design and purpose</i>
	7.2 <i>Special features</i>
	7.3 <i>lugano in parts</i>
	7.4 <i>The lugano ready for use</i>
	7.5 <i>Assembly of Electrical Braking System (EBS)</i>
	7.6 <i>Control</i>
	7.7 <i>Disassembly</i>
	7.8 <i>Relocation</i>
	7.9 <i>Functioning advice</i>
	7.10 <i>Disposal</i>
	7.11 <i>Troubleshooting</i>

