

livorno low

The model for comprehensive well-being



Numerous equipment variants as well as a large selection of colours and end panels make the **livorno low** a **flexible nursing care bed** for rehabilitation and care facilities.

In addition to its elegant design, it is able to convince with practical functions and high suitability for everyday use.

The livorno low cuts a good figure in every respect: **economically** thanks to its attractive price-performance ratio and **practically** due to its functionality, both to a lasting relief in the daily care.

Last but not least, a **robust construction and durable materials** ensure that the livorno low is an investment that will surely pay off.



Mona (35), Passau

“The livorno low is a real all-rounder, because it can be optimally adjusted to every resident and offers all the functions you need in everyday life”.

livorno
LOW

livorno
LOW

livorno low



Suitable for TEXIBLE Wisbi® Page 42

| Dimensions | |
|-------------------|------------------|
| Outer dimensions | 105 x 212 cm |
| Lying surface | 90 x 200 cm |
| Special widths | 100, 110 cm |
| Special lengths | 190, 210, 220 cm |
| Height adjustment | 29 - 73 cm |
| Lifter clearance | 15 cm |

| Weights | |
|--------------------|------------|
| Bed weight approx. | ca. 109 kg |
| Safe working load | 200 kg |
| Max. person weight | 165 kg |

| Side rails | |
|----------------|---|
| Material: Wood | |
| Version | continuous long divided short divided head-section |

| | |
|---------------------|---------|
| Height*: | |
| - With wooden slats | 39.5 cm |
| - With ripolux® neo | 35 cm |

Further equipment features

| | |
|------------|---------------------|
| End panels | I, IV, V, VII, VIII |
|------------|---------------------|

Motor
 Splash-proof according to IPX4, high-quality spiral cable with strain relief and kink protection, SMPS plug transformer with extra-low safety voltage offering safety from the socket

*Distance between upper edge of lying surface and upper edge of the raised side rail



Product specifications, end panels and accessories can be found from page 56 onwards



Hand control with locking function



Central braking, optionally with aluminium covered castors



Mobilisation-supporting side rails

The product advantages of livorno low:



End panel Typ I

- ⊕ Central braking
- ⊕ Lifter clearance
- ⊕ Low bed entry height 29 cm
- ⊕ Hand control with locking function
- ⊕ Various comfort lying surfaces
- ⊕ Various side rail options

More information:

www.bock.net/x39a-en

End panel type V with curved handle

Wooden slats with tension regulators

4-section lying surface

Hand control with locking function

Aluminium covered castors

Vertical height adjustment

Short divided side rails

Central braking