


practico 25|80

Anything but normal ...



Practical functions, a generous **height adjustability**, **comfort and safety** as well as an appealing design, which can be perfectly integrated into many living environments: Even the **practico 25|80** looks like a real “bock” at first glance - a fully developed nursing care bed with the ingenuity of many decades of experience and inspired by the people who work in the care sector day after day. In addition to the **bed extension** (side rail exchange required), it has many other useful features, such as central braking,

head and leg down position, or the **comfort sitting position**.

In addition, the practico 25|80 can be equipped with attractive extras, ranging from the elegant aluminium covered castors, to underbed lighting and various wall distance holders. For an aluminium version with integrated, tool-free bed and side rail extension please have a look at the previous pages where the practico alu plus and the practico alu 25/80 are described.



Eileen (42), Bad Lippspringe

“Everything included: The practico 25|80 has everything a nursing care bed needs - and it’s also super comfortable”.

practico
25|80

practico
25|80



Suitable for TEXIBLE Wisbi® Page 42

Dimensions	
Outer dimensions	105 x 210 cm
Lying surface	90 x 200 cm
Special widths	100, 110 cm
Special lengths	210, 220 cm
Height adjustment	25 - 80 cm
Lifter clearance	15 cm

Weights	
Bed weight approx.	147 kg
Safe working load	220 kg
Max. person weight	185 kg

Side rail	
Material: Wood	
Version	continuous long divided short divided head-section
Height*:	
- With wooden slats	39,5 cm
- With ripolux® neo	35 cm

Further equipment features	
End panels	I, IV, V, VII, VIII
Motor	
Splash-proof according to IPX4, high-quality spiral cable with strain relief and kink protection, SMPS plug transformer with extra-low safety voltage offering safety from the socket.	

*Distance between upper edge of lying surface and upper edge of the raised side rail



Product specifications, end panels and accessories can be found from page 56 onwards




Height adjustment to the optimum care height of 80 cm



Long divided side rails



Central braking of all four castors optionally with aluminium covered castors



Hand control with buttons for a comfort sitting position and bed lighting



More information:



www.bock.net/x37k-en

The product advantages of practico 25|80:

- ⊕ Height adjustability from 25 - 80 cm
- ⊕ Central braking with directional giving device
- ⊕ Lifter clearance 15 cm
- ⊕ Comfort sitting position
- ⊕ Different side rail variants to choose from
- ⊕ Integrated bed extension (side rails must be exchanged)

