

livorno low

The elegantly designed, multifaceted nursing care bed offers a range of customised configuration options

Perfectly suited for care.

.bock[®]///





livorno low

A wide range of equipment

The attractive livorno low nursing care bed is available with a **wide range of equipment**, and its elegant design and user-friendly functionality makes it easy to switch to a new bed.

Even the basic version of the livorno low includes everything a good nursing care bed needs, from an adjustable height to space for lifters.

The **electric adjustment functions** are particularly easy to use, enabling you to easily select a comfortable height for getting into or out of the bed or a pleasant resting position using the hand control.

The livorno low can be brought down to a particularly **low level for getting into the bed**, which also enables the resident to sit safely and comfortably on the edge of the bed.

Attractive equipment options such as covering of castors and a large selection of colours makes the livorno low a practical and comfortable bed, and one that is highly regarded in many care homes and rehabilitation centres.

Divided side rails on the head end offer additional stability and safety for the resident when getting out of bed. Under-bed lighting ensures that residents can safely leave and return to the bed at night.



Double-sided continuous side rails



Long-divided side rails



For more information on the livorno low, visit
www.bock.net/livorno

Configuration options

You can choose from different colours for the livorno low. Practical extras allow you to perfectly adapt the livorno low to your needs.

What the livorno low offers:


- /// Height-adjustable from 29 to 73 cm
- /// Central braking
- /// Lifter space > 15 cm
- /// Hand control with locking function
- /// Wide range of comfortable lying surfaces
- /// Variable side rail solutions
- /// Approved for use as a medical product



The optional aluminium castor covering makes the livorno low especially comfortable.




The hand control has large buttons and is easy to use.



The optional under-bed light can be operated using the hand control, just like the reading lamp.




Your choice of colours:



Beech, steamed R24034



Havanna D7935




Sonoma oak, light R20128




Beech R24015



Cherry R42009



Maple, dark R27017



Walnut, dark R30011

Your choice of end panels:



End panel type I




End panel type V



End panel type VII

Specifications

External dimensions:	105 x 212 cm
Safe working load:	200 kg
Max. person weight:	165 kg
Height adjustment:	29 – 73 cm
Lifter space:	15 cm
Side rail height:	40 cm with wooden slats / 35 cm with ripolux® neo
Lying surface dimensions:	90 x 200 cm
Lying surface:	4-section



Electric height adjustment and electric back and leg section adjustment with automatic triple-action function

Top quality made in Germany

All bock beds are approved for use as medical products and are designed and manufactured at our location in Verl, Germany.

Our beds are available worldwide at specialised medical retailers. You can find your local retailer at <http://en.bock.net/retailer>.



Further information on bock's comfort and nursing care beds is available on our website. Simply scan this QR code and you will automatically be brought to our website.

.bock

.bock[®] ///

Hermann Bock GmbH
Nickelstraße 12
33415 Verl
Germany

Phone: +49 (0)5246 9205-0
Fax: +49 (0)5246 9205-25
E-mail: info@bock.net
Website: www.bock.net

© registered trademark
Order number: 890.02396
As of October 2018
Subject to technical changes.