

The latest generation of beds for care at home.



The new bock domiflex® 3 takes an important step into the future to improve care at home. For retailers, caregivers as well as for people in need of care the domiflex® 3 is the ideal bed for care at home. Even the standard model has an extended height adjustment range of 35–80 cm, and offers an increased patient weight of up to 155 kg and thus a higher safe working load of 190 kg.

domiflex® 3	
Dimensions	
External dimensions	103 × 210.7 cm
Lying surface	90 × 200 cm
Special dimensions	Bed extension kit for length 220 cm available
Height adjustment	35–80 cm
Lifter space	> 15 cm
Weights	
Max. patient weight	155 kg
Safe working load	190 kg
Bed weight	78 kg
Max. weight per component	< 18 kg
Motoren	
Splash-proof according to IPX4, high-quality spiral cable with strain relief and kink protection, SMPS plug transformer with extra-low safety voltage offering safety from the socket.	



You can find all the advantages of bock domiflex® 3 here:

Lightweight components and a tool-free connection technology ensure a quick assembly and disassembly

4-section lying surface with comfortable wooden slats

35 cm for a safe entry and exit

domiflex[®] 3

All advantages at a glance



Height adjustment range 35–80 cm

- ⊕ Extended height adjustment range of 35–80 cm
- ⊕ Ergonomic and comfortable entry and exit at 35 cm
- ⊕ Back-friendly care height, adjustable up to 80 cm



Choice of two colours

- ⊕ Classic: Steamed Beech
- ⊕ Modern: Sonoma Oak Light



Cable-connected hand control

- ⊕ With locking function and numerous adjustment possibilities
- ⊕ Soft buttons provide comfortable haptics
- ⊕ Safe handling



Wide range of accessories

- ⊕ More comfort, for instance thanks to side rail bumpers
- ⊕ More safety, for instance thanks to the grab rail that offers an easy mobilisation



Transport and storage system

- ⊕ The space-saving transport and storage system as well as lightweight components enable a simple handling



Tool-free connection

- ⊕ The tool-free connection technology enables a quick assembly and disassembly



Hermann Bock GmbH
Nickelstrasse 12
33415 Verl
Germany

Phone: +49 5246 9205 0
Fax: +49 5246 9205 25
Internet: www.bock.net
Email: info@bock.net

.bock[®]